




It is our job to make
your glasses
comfortable

Since 1977

Hasegawa Bicoh co., ltd.



Thank you very much for your continued patronage of our products.
We, Hasegawa Bico Co., Ltd., established in July 1977 with our policy
"It is our job to make your glasses comfortable".
We have developed and manufactured
"eyeglass related products not found anywhere in the world".

Thanks to everyone, we were able to celebrate our 40th anniversary in July 2017.
Again, I am sorry for the gratitude and support of everyone
who had patronage for many years.
We are good at silicone products, and products such as "Megalock" and "Cell Seal"
The world's first "Silicone nose pad", "Cushion nose pad",
"Silicone temple tip" and other glasses parts,
Especially for the trial frame, it occupied about 80% of domestic market share.
We are also exporting to more than 20 overseas countries.
The society of the future will become international era,
the era will be more demanding speed than ever before.
Among them, we respond "seriously" and "promptly"
to opinions and requests from customers daily,
We are committed to working on the development of better products.
In 2020, the pandemic of the COVID-19 faced an unprecedented crisis
in the world and had a significant impact on the eyewear industry.
However, we will continue to increase our potential and
move forward to develop new products every day,
so we appreciate your guidance and encouragement.

President Isao Kudo

CONTENTS



- 2 CELL SEAL U
CELL SEAL M

- 3 SPONGE PUFF NOSE PAD
UP & DOWN NOSE PAD



- 4 MEGALOCK V
MEGALOCK F
MEGALOCK M



- 5 MEGANE GRIP
ROUGLASS L



- 6 CLEAN COAT



- 7 BiBa



- 8 SILICONE NOSE PAD DX
CUSHION NOSE PAD C-one
CUSHION HILON NOSE PAD



- 9 HILON NOSE PAD EX
HILON NOSE PAD S
HILON NOSE PAD R
CUSHION NOSE PAD AR



- 10 DRIVER(+) FOR REMOVING
BROKEN SCREWS

- 11 GACHIDORA



- 12 CSL PLIER



- 13 CS PLIER
CELL MORI
SILICONE COAT SET



- 14 SILICONE TEMPLE TIP (ROUND HOLE)
SILICONE TEMPLE TIP (FLAT HOLE)
SILICONE TEMPLE TIP (SLIM TYPE)

- 15 SOFURE TEMPLE TIP
SILICONE CURL TIP FOR ADULT
SILICONE CURL TIP FOR KIDS



- 16 UNIVERSAL TRIAL FRAME
580Univ Pro
580Univ Pro – LO

- 17 580Univ
580Univ – LO

- 18 MC674
MC486

- 19 FIX TYPE TRIAL FRAME
SIMPLE BC
SIMPLE BC – LO
COLOR SIMPLE BC



- 20 TRIAL SET FOR HIGH CURVE




- 21 FITTING MEASURE
CLUTCH GLASS α
for "EYELID PTOSIS"

All made in Japan



To prevent slip down of the plastic frame

CELL SEAL® U



- Just paste on the nose pad.
- When the nose pad is raised, the distance between the eyes and the lens is separated. Eyelashes and cheeks are less likely to touch the eyeglasses.
- The shape of the cell seal has been changed to make it easier to turn up. Easy to peel off and paste by making cuts in the release paper.
- It is transparent color, made of silicon which is safe to the skin.



PAT no. 3157348

Size (Thickness)	S (1.3mm) / M (1.8mm) / L (2.5mm) / LL (3mm)
Material	Silicone
Package Box	1 pair Including 10 pairs.

How to use

Raise the nose pad




Paste the nose pad of the frame sandwiching.

Narrow the nose width




Cut and paste the reinforcement part.



Small Size CELL SEAL

CELL SEAL® M



- For small plastic frame.
- Easy to turn up and paste.

PAT no. 3157348

Size (Thickness)	M (1.8mm) / L (2.5mm) / LL (3mm)
Material	Silicone
Package Box	1 pair Including 10 pairs




It is very fluffy

PAT no. 3202187



SPONGE PUFFS NOSE PAD®

- The shape of the round dome
- Antibacterial sponge puff
- Use a special medical adhesive seal
- Strong adhesive power



SALES POINT

Keep the make up of nose
Glasses do not feel heavy.



Just paste on the nose pad of glasses.



Others



Ours Thickness 2mm

Package	2 pairs
Box	including 10 package
Material	Nitrile Butadiene Rubber




PAT no. 4833954




For the progressive lens use

UP & DOWN NOSE PAD®



- The nose pad moves up and down by 5 mm. It makes it easy to see nearby.
- V character shape. Vertically up and down. The nose width does not change.
- Silicone surface



Since the near portion of the lens approaches the eye, You can see small characters and others nearby.

Package	1 pair
Type	Screw in / Push on
Size	15.5 × 8mm

Display set	
5 pair set	12 pair set
Screw in	Screw : 10 / Push : 2
Hole size	1.2 φ
Material	Nylon



Suitable for various size temple tip

MEGALOCK® V



- Just insert to the temple tip. Glasses are fixed.
- Especially recommended for those who like to move such as sports.
- Since the distance between the eye and the lens is fixed, the burden on the eye is reduced.



The adjust parts for the slim temple tip.


* Width of temples : about 5mm - 10mm

Package	1 pair	Include adjust parts & case
Color	BLACK / BLUE / GREEN / PINK	
Material	Silicone	
Box :	1) 1 color : 10 pair 2) Assortment : 10 pair (Black:3 / Blue:3 / Green:2 / Pink:2)	



Large size Megalock

MEGALOCK® F




* Width of temples : about 4mm - 8mm

Package	1 pair	Include adjust parts & case
Color	BLACK / BROWN / BEIGE	
Box	10 pairs	
Material	Silicone	



Small size Megalock

MEGALOCK® M



* Width of temples : about 4mm - 8mm

Package	1 pair	Include adjust parts & case
Color	BLACK / BROWN / GREEN	
Box	10 pairs	
Material	Silicone	





Just cover the temple tip by the silicone tube.

MEGANE GRIP



- Just insert to the temple tip.
- Suitable for temple tip because stretch silicone.
- Safety silicone material.



- (M) Metal frame
- (L) Plastic frame




Package	1 pair	Box	10 pairs
Color	Brown	Material	Silicone
Size	M / L	Outer/inner	M : 2.6 φ / 4.3 φ L : 3.1 φ / 4.8 φ


Change to bifocals just by attaching to glasses.

ROUGLASS® L

The licence of "manufacture and sales of medical equipment"
No. 13B3X00177



- It can be installed with one touch.



Just pinch the clip and insert it inside the eyeglasses. Add the lens power to make it easy to see nearby.

Sales point

- I can not see detailed characters clearly.
- For nearness a little bit hard to see, but a bit more.
- I want to use my current glasses.
- It is troublesome to carry reading glasses on a trip or going out.

Reference for lens power

AGE	POWER
40-50	+1.00
45-55	+1.50
50-60	+2.00
55-65	+2.50
60-70	+3.00
70 ~	+3.50



Power	+1.00 / +1.50 / +2.00 / +2.50 / +3.00 / +3.50
Material	Lens : Acrylic Clip : Nylon
Display	9 pieces(+1.00 :1 / +1.50 :2 / +2.00 :2 / +2.50 :2 / +3.00 :1 / +3.50 :1)



Fluorine blending in lenses, slippery, long-lasting glasses

CLEAN COAT

Make a coating film of fluorine and protect glasses from scratches and dirt.

Also for smart phones, PC monitors and other uses other than glasses.



Bulk	10mℓ
Material	Surfactant · Fluorine oil

Both cleaners and coatings can be done!
Alcohol is not used at all.

Effect on lens

- Removes a wide variety of dirt on the entire lens surface.
- By making a coating film of fluorine on the whole lens surface
 - Prevent contamination and scratches.
 - Easily remove dirt.
 - Touch feeling when wiping the lens refreshing It is slippery.

Effect on frame

- By making a coating film of fluorine on the whole frame
 - It protects from hair tonic for Temple and temple end.
 - If you are allergic to the rash in Temple Prevent.
 - Prevention of whitening of plastic and tortoise shells.



Package	For sale : 10mℓ (20 pcs.)	For store demonstration : 200mℓ (3 pcs.)	Replenishment : 1ℓ (1 pcs.)
---------	------------------------------	---	--------------------------------









Optical frame for a baby (1-2 age)

BiBa

BiBa

Emphasis on safety · cuteness, hold tolerance for young children.
Since after parts are also available, various customization is possible.
Please try the shape that you can feel comfortably without disliking glasses.

	
①+②+③+⑦	①+②+④ / ⑤+⑥
Head band type	Curl tip type
	
①+②+④ / ⑤+⑥+⑦	①+②+④ / ⑤+⑦×2本
Curl tip & Head band type	T-style type
	
①+②+④ / ⑤+⑦×2本	①+②+④ / ⑤+⑦×3本
Goggle type	X-style type











Curl tip type
(Blue)



Head band type
(Pink)

Size	Frame : 36 / 38
Color	Pink / Blue

After part

① Front frame	② Silicone nose pad	③ temple for the band	Parts set
			
④ short temple ⑤ Long temple	⑥ Curl tip	⑦ Band	
		Front  Back 	


Color : Pink / Blue

All made in Japan



Standard of the silicone nose pad

SILICONE NOSE PAD DX




Size	L (20mm) /M (17mm) /S (13mm)		
Color	Pad : Transparent	Hole size	1.2 ϕ
	Metal: Gold or Silver		
Package	1Box : 9 pairs (Gold :6 pairs/Silver :3 pairs) (each size)		
	Bulk : 50 pairs (each size)		

Gently absorb the movement of eyeglasses by the cushion function

CUSHION NOSE PAD C-one

PAT. no.3208475


- A resin cushion mechanism is provided in the conventional core portion, By integrally molding (making all shapes with one material) It enables unprecedented strength, springiness and flexibility We realized comfortable wearing feeling.
- Newly developed material Using ethylene elastomer, it is excellent in durability and chemical resistance, It is also used for medical use (clear food sanitation law), it is also a gentle pad on the skin.



Size	L (17mm) /M (14.6mm) /S (11.6mm)		
Color	Transparent	hole size	1.1 φ
Type	Screw in / Push on		
Package	1 Box :25 pairs (sizes and types each)		
	Bulk :50 Pairs (sizes and types each)		

Clear cushion pad of special plastic

CUSHION HILON NOSE PAD



Size	L (17mm) /M (14mm)		
Color	Transparent	Hole size	1.1 φ
Package	1 Box :9 pairs (each size)		
	Bulk :50 pairs (each size)		



- High performance nose pad using special engineering plastic.
- Nose pad traces are difficult to attach to the nose.
- Cosmetics and sebum dirt are difficult to adhere.


Hilon nose pad

SERIES

There is cushioning property by adoption of silicon ring

HILON NOSE PAD EX

PAT. no. 3169961



- Soft cushioning.
- No sound and no clattering.

Size	M (14.4mm)/S (13mm)
Hole size	1.2×1.6mm (1.2 screw available)
Color	Pad :Clear beige / Metal :Gold
Package	1 box :9 pairs (each size) Bulk :50 pairs (each size)



Standard Hilon pad

HILON NOSE PAD S

Size	14.4mm	Hole size	1.2φ
Color	Pad :Clear beige / Metal :Gold		
Package	1 box :9 pairs / Bulk :50 pairs		



Hilon pad for Men

HILON NOSE PAD R

Size	15mm	Hole size	1.2φ
Color	Pad :Clear blue / Metal :Silver		
Package	1 box :9 pairs / Bulk :50 pairs		

Cushion function of the new mechanism realized

CUSHION NOSE PAD AR

PAT no. 4833790



This gap is cushion

View from top side.

- Adopted comfortable material "ethylene elastomer".

Size	14mm		
Color	Transparent	Hole size	1.1 ϕ
Package	1 box : 9 pairs / Bulk : 50 pairs		



Easy to remove broken screws

Driver(+) for removing broken screws

NEJI TORU[®]plus ⊕ PAT.P

Bicoh's
Eyeglass tools

- A new screwdriver that can remove broken screws.
- This driver is dedicated screwdriver to remove with a crushed screw head.
- For hinge screw, rim screw, pad screw.
Phillips(+) screw only. (Depends on size and condition.)




Actual size


length :122mm
weight:46g

Try this before sending for repair!

Can remove a crushed screw.



- ⚠ This driver is for counterclockwise turn only.
- ⚠ Do not turn clockwise.



• Hinge screw



• Pad screw



• Rim screw

Screw diameter	1.0mm 1.2mm 1.4mm
Package	Driver 1 pc. Wrench for replace the driver tip.

Optional replacement parts

bit for the driver



Actual size



- A screwdriver that fits very well in the groove of the screw.
- It doesn't move rattlingly.
- The power to turn is efficiently transmitted to the screw.

Attention




Silicone tube for
tip protection to pick up
screws.

Attention

Precision Phillips(+) screwdriver

The tip of the screwdriver
was designed only for eyeglass screws.



Actual size

length :115mm
weight:33g



• Hinge screw



• Pad screw



• Rim screw

Optional replacement parts

bit for the driver



Actual size

Screw diameter	1.0mm	1.2mm	1.4mm
Package	Driver 1 pc. Wrench for replace the driver tip.		



The plier can adjust various nose pad arms.

CSL PLIER

PAT no. 3157348

Bicoh's
Eyeglass tools

2020
"The most emotional eyewear"
Grandprix




Attention


- CSL plier can adjust "click-in type" nose pad arm smoothly.
- Good hold the metal in the nose pad part with "grooves inside the tool" it can be adjusted with small power.
- It can also be used for "screw-type" pad box.

New Function of CSL plier

(CS plier do not have)




Very useful for adjusting these types of Clings.




New Function of CSL plier

(CS plier do not have)



Vertical groove at the tip.

- It also has the function of a CS plier.



Please grasp the screw so that the screw is in the hole of the plier tip.

Accessory parts



A (fix for end piece)



B (fix for temple)




Insert to the top of the plier.

Package

Plier 1 pc.
Accessory A (fix for end tip) B (fix for temple)



Bicoh's Eyeglass tools



The plier for adjust of nose pad

CS PLIER

- The CS plier is useful tool. Adjust nose pad box easier.
- Since you can grasp only the pad box surely, There is nothing to hurt.
- You can also use it for regular pads. (silicon nose pad, hard nose pad, etc.)
- Include accessory parts. (same as the CSL Plier.)

Package

Plier 1 pc.
Accessory A (fix for end tip) B (fix for temple)

For celluloid and acetate



No.1

No.2

No.3

- Please bond with "Silicone Coat Set".
- 1) Take off the oil stains on the part to be coated.
 - 2) When partially painting in a place without a step please mask with a cellophane tape etc.
 - 3) Coat the primer and let it dry for more than 5 minutes.
 - 4) Have the coat part level and apply silicone with a brush.
 - 5) Paste the cell mori, touch lightly with your fingers and pull out the air.
 - 6) Dry about 2-10 hours.

Silicone nose pad for the plastic frame

CELL MORI



Size	No.1 (1mm) / No.2 (1.5mm) / No.3 (2.2mm)
Package	1 box : 12 pairs (per size)
Material	Silicone



Primer



Silicone glue



Coating brush

It can adhesiving and coating

SILICONE COAT SET

- Apply to the nose pad or temple tip part, etc. to prevent pain and gaps.

Functions

- Protect the human skin from plating, apply to the parts to prevent rash.
- Replacement of silicone portion of silicone processed frame, for repair bonding.
- It can also be used for coat for nose pad of plastic frames.


For celluloid and acetate

Package

Silicone glue×1 pc. / Coating brush×1 pc. / Primer×1 bottle



SILICONE TEMPLE TIP (ROUND HOLE)



Color	Black /Brown
Hole size	1.4 ϕ — 1.5 ϕ (round hole)
Package	1) 5 pairs (Black :2 pairs / Brown :3 pairs) 2) 1 color 5 pairs (per color)
Material	Cellulose propionate, Silicone

SILICONE TEMPLE TIP (FLAT HOLE)



Color	Black / Brown / Demi
Hole size	1.4 flat — 1.5 flat
Package	1) 5 pairs (Brown:2 / Demi:2 / Black:1) 2) 1 color 5 pairs (per color)
Material	Cellulose propionate, Silicone

SILICONE TEMPLE TIP (SLIM TYPE)




Color	Black / Brown / Demi / Pearl beige / Transparent
Hole size	1.2 ϕ — 1.3 ϕ (round hole)
Package	1) 5 pairs (per color) 2) 1 color 5 pairs (per color)
Material	Cellulose propionate, Silicone



Soft material temple tip for round and flat hole

SOFURE TEMPLE TIP



- Polyester elastomer (soft plastic).
- Very durable (cracking is also difficult to occur)
It is gentle for the skin.

Color	Black / Brown
Size	1.3φ — 1.4φ (round and flat)
Package	1) 8 pairs (Black :3 pairs/Brown :5 pairs) 2) 1 color 8 pairs (per color)
Material	Polyester elastomer

SILICONE CURL TIP FOR ADULT




SILICONE CURL TIP FOR KIDS



- Soft spring properties with special resin.
- Special silicone on the surface gentle to your skin and soften your ear hit.
- It can be adjusted according to the shape of the ear with hot air.
- Abundant color variation.

Color



Color	Black / Brown / Dark red / Red / Green / Dark blue / Violet/ Light blue	
Size	1.3φ — 1.5φ round hole	
Package	for Adult : 1) 4 pairs (Black :1/ Brown :2 / Dark red :1)	2) 1 color 4 pairs (per color)
	for Kids : 1) 4 pairs (Light blue / Green / Red / Violet :each 1 pair)	2) 1 color 4 pairs (per color)
Material	Silicone, Polyester elastomer	



Add new function to 580 Univ. Lightest 53g ※1

580 UniV® Pro

PD 50 – 80

PAT.P.

The licence of “manufacture and sales of medical equipment” 13B3X00177

Attention

- Nose pad moves front and back by 5mm
- With corneal vertex distance measurement scale
- 5 Lenses (Front:3pcs. & back:2pcs.)

2019

“The most emotional eyewear”
Grandprix



Half rim type for progressive lens use

580 UniV® Pro-LO

PAT.P.

The licence of “manufacture and sales of medical equipment” 13B3X00177

- Half Rim type
- Lightest 51g ※1




- PD 50~80 mm (25-40 mm per right and left)
- Adjust : Temple length · Height & width & open angle of nose pad · Front angle




PD	50 – 80mm (25 – 40mm per right and left)		
Nose pad	Height : 9mm (Stepless), width & open angle by 2 actions, Front and back : 5mm (Stepless)		
Temple	Length : 30mm (per 5mm) Front angle : up & down 7 degrees		
Astigmatic axis	Rotation function (5 degree)	Lens	Front : 3 pcs. / Back : 2 pc.
Weight	53g	Color	Grey

※1 Compared with the trial frame (left and right PD move independently, Nose pad moves front and back) by other companies.


New Functions



Turn the knob on the front
Nose Pad moves front and back by 5mm.



With corneal vertex distance measurement scale



5 Lenses (Front:3 pcs. & back:2 pcs.)



The left and right PD move independently. Lightest 52g^{※2}

580 UniV[®]

PD 50 - 80

PAT.P.

The licence of "manufacture and sales of medical equipment" 13B3X00177



- PD 50~80 mm (25-40 mm per right and left)
- Adjust : Temple length • Height & width & open angle of nose pad • Front angle
- 4 Lenses (Front:3 pcs. & back:1 pc.)


Half rim type for progressive lens use

580 UniV[®]-LO

PAT.P.

The licence of "manufacture and sales of medical equipment" 13B3X00177







- Half Rim type
- Lightest 49g^{※2}



PD	50 - 80mm (25 - 40mm per right and left)		
Nose pad	Height : 9mm (Stepless), width & open angle by 2 actions		
Temple	Length : 30mm (per 5mm) Front angle : up & down 7 degrees		
Astigmatic axis	Rotation function (5 degree)	Lens	Front : 3 pcs. / Back : 1 pc.
Weight	52g	Color	Grey

※2 Compared with the trial frame (left and right PD move independently, Nose pad moves front and back) by other companies.

Instruction






universal trial frame with one-touch PD adjustment

PAT no. 4121660

MC 674

PD 54 - 74

The licence of "manufacture and sales of medical equipment" 13B3X00177



PD	54 — 74mm		
Nose pad	Height, width, open angle by 2 actions		
Temple	Length : 30mm (per 5mm) Front angle : up & down 7 degrees		
Astigmatic axis	Rotation function (5 degree)	Lens	Front : 3 pcs. / Back : 1 pc.
Weight	50g	Color	Grey

Full rim & Straight temple

674



Half rim & Curl temple

674



Half rim & Straight temple

674




Universal trial frame for women and kids

PAT no. 4121660

MC 486

PD 48 - 60

The licence of "manufacture and sales of medical equipment" 13B3X00177



PD	48 — 60mm		
Nose pad	Height, width, open angle by 2 actions		
Temple	Length : 20mm (per 5mm) Front angle : up & down 7 degrees		
Astigmatic axis	Rotation function (5 degree)	Lens	Front : 3 pcs. / Back : 1 pc.
Weight	48g	Color	Grey

Full rim & Straight temple

486



Half rim & Curl temple

486




Half rim & Curl temple

486



Spare for temple tip



• Straight temple tip



• Curl temple tip

Package	1 pair
For	674 series, 486 series 580 series



The licence of "manufacture and sales of medical equipment"
13B3X00177

Standars PD fix type TRIAL FRAME

SIMPLE BC

PD	54	56	58	60	62	64	66	68	70	72	74
----	----	----	----	----	----	----	----	----	----	----	----



- Ultra lightweight (35 g) with few breakdown frames.
- Movable temple tip and gently embracing your head.
- Safe hinge structure without risk during desorption.

Simple BC for KIDS

COLOR SIMPLE BC

PD	48	50	52
----	----	----	----



* Discontinued manufacturing : Pink (Color simple, spare temple tip)

PD Size	48 / 50 / 52
Temple	Length: 9mm (per 3mm) *PD52 : 12mm (per 3mm)
Astigmatic axis	Rotation function (5 degree)
Lenses	Front: 3 pcs.
Weight	33g
Color	Blue


For the proffessive lens use

SIMPLE BC-LO

PD	56	58	60	62	64	66	68	70	72
----	----	----	----	----	----	----	----	----	----



PD Size	54 (simpleBC) / 56/58/60/62/64/66/68/70/72/74 (simpleBC)
Nose Pad	Hight 2 action
Temple	Length 12mm (per 3mm)
Astigmatic axis	Rotation function (5 degree)
Weight	35g
Lenses	Front: 3 pcs./Back: 1 pc.
Color	Black



BLUE



BLACK


Spare temple tip

for Simple series

other option parts for our trial frame

Clutch function for the trial frame

For examination of ptosis



Head band

Prevents slip of the trial frame

* Some older trial frames,
It may not be installed.





Feel free to experience how the "high curve lens" looks


TRIAL SET FOR HIGH CURVE

PAT no. 3178090

- For customers considering sunglasses with high curve
"Prescription glasses"



Just by attaching a high curve lens to the frame
You can experience how the high curve looks.



Lens astigmatism inside
Spherical power lens outside



- It is adjustable for customer PD.
- Put an astigmatic lens inside.
- * Since high-curve lenses are not included, please arrange at your shop.

Package

A) Frame	: 1 pc	B) Seal for lens	: 1 pair
C) Ring for astigmatic	: 1 pair	D) Scale for axis degree decision	: 1 pc
E) Ring for high curve	: 3 pairs / 6 pairs / 10 pairs		



Tool for measuring head size

FITTING MEASURE

PAT no. 3191368

- For frame pre-fitting
- Temple opening width and temple tip bending length are known by mm.
- It is similar to eyeglasses.
- You can appeal your skills.



Package 1 pc.



We are accepting OEM for trial lens set.



For those who need help with ptosis


NEW

CLUTCH GLASS α for “EYELID PTOSIS”

alpha

- Clutch glass alpha, with elastic spring, with special plastic, soft eyelid fit.
- Eliminate the discomfort of vision disorder caused by ptosis of the eyelid.
- It can be attached to a metal frame with a temple width of 5 mm or less.

△ Use “CLUTCH GLASS II” for plastic frame.



Point for improvement

Durability has improved, because by adding a one-turn twist at the base.

Package Sold separately right and left

All made in Japan

Slip off of the eyeglasses causes fatigue and eye strain.



2021.10

 *Hasegawa BicoH co., ltd.*

2-41-4 Matsumoto, Edogawa-ku Tokyo 133-0043 Japan

Tel : 81-3-3655-4321 Fax : 81-3-3655-4472

<http://www.bicoh.co.jp>

E-mail : bicoh-export@bicoh.co.jp